

PRODUCT LIST FOR TREATMENT OF VITILIGO / LEUCODERMA AND PSORIASIS

Features:

- Equipped with UVA/NBUVB tubes as per doctors suggestion.
- Compact & Portable.
- Used to treat localized area of the body.
- Less space required.
- Ideal for clinics.

Home Therapy Unit:

- Option A: 2 feet UVA tubes (4 Nos) with timer.
- Option B: 2 feet NBUVB tubes (4 Nos) with timer.



Single Panel UV Therapy Unit:

- Option A: 6 feet UVA tubes (6 Nos) with Timer & Dosimeter for setting jules.
- Option B: 6 feet NBUVB tubes (6 Nos) with Timer & Dosimeter for setting jules.



Hand & Foot Unit:

- Home treatment of both Hand and Foot simultaneously.
- 4 panels of 2 feet tubes (8 nos UVA + 8 nos NBUVB).
- Timer & Dosimeter for setting jules.



Spot PUVA:

- With shutter to focus on small spots.
- UVA/NBUVB tube (6 inch - 1 no).



Portable PUVA:

- Handheld unit .
(6 inch UVA tubes - 2 nos.) .



Hand NBUVB:

- Without Timer (NBUVB):
- Comb attachment for scalp treatment.



Precaution: Use UVA protective goggles (Do not see tubes without goggles)
Dosage: as per Doctors recommendation.

Authorised Distributors :

JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072, e-mail : jollytechnocrats@gmail.com, jollytechnocratsservice@gmail.com
website : www.jollytechnocrats.in, Mobile : +91-9650595508

