

Fluid Warming Cabinet

JOLLY TECHNOCRATS

An ISO 9001:2015 Company



Model : IV Care 45



Model : IV Care 100

TECHNICAL SPECIFICATION

- ❖ Bright LED Display with 0.1' C resolution
- ❖ Microprocessor based digital monitor to display :- Temperature Reading Shown on Double Display.
- ❖ Temp. Reading Shown in °C / °F , Temp, value setable $\pm 25^{\circ}$ to 80°
- ❖ Audio / Visual Alarms for :- Temperature, door open, power failure, sensor failure
- ❖ LED Indications with :- temp, fail power fail, door open, Heat on, Silence
- ❖ Integrated Voltage protective device against excessive voltage fluctuation.
- ❖ Double vaccumed glass door as standard, solid door (optional)
- ❖ Fully St. Steel body & size are also available as per customer requirement.
- ❖ Password Protection to enter Main Setting
- ❖ Inner chamber of Stainless Steel for rust free interior
- ❖ Outer Chamber choice of **A.** Powder Coated Material, **B.** Stainless Steel 202 Grade
C. Stainless Steel 304 Grade.

Index for Fluid Warming cabinet			
Model	Inner Capacity (Approx)	Storage Capacity	Size Approx
I V Care - 45	145 Lts	10-15 Bottle of 3 Liter	5 Cubic Feet
		40-50 Bottle of 1 Liter	
IV Care - 75	205 Lts	25-30 Bottle of 3 Liter	7 Cubic Feet
		70-80 Bottle of 1 Liter	
IV Care - 100	285 Lts	40-50 Bottle of 3 Liter	10 Cubic Feet
		100-110 Bottle of 1 Liter	
IV Care - 125	345 Lts	60-70 Bottle of 3 Liter	12 Cubic Feet
		120-130 Bottle of 1 Liter	



JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072, e-mail : jollytechnocrats@gmail.com, jollytechnocratsservice@gmail.com

website : www.jollytechnocrats.in, Mobile : +91-9810025581, +91-9650595509