

# Looks 4000 A / 5500 A / 7500 A

**NEW**



“No head bent loupes”



## ❖ DESCRIPTION

- This product has the characteristic that reduces the stress applied to the cervical and lumbar spine when performing medical treatment by bending the head less than the conventional product.
- The telescope of KEPLERIAN TYPE which PRISM is contained internally offers a height resolution.
- This is product to make efforts to overcome the defect to be occurred in high magnification as follows - FOV to narrow, sharpness to fall, and weight to increase etc.
- We recommend to the person with the experience which used lower magnification.

## ❖ FEATURE

Model	Looks 4000A	Looks5500A	Looks 7500A
Total weight	55.1 g	59 g	68.1 g
Field of view (FOV)	100 mm @500 mm	73 mm @500 mm	53 mm @500 mm
Resolution	More than 28.5 lp/mm	More than 28.5 lp/mm	More than 28.5 lp/mm
Image distortion	Less than 1 %	Less than 1 %	Less than 1 %
Image flatness	Less than 0.2 diopter	Less than 0.2 diopter	Less than 0.2 diopter