



Looks® SURGICAL TELESCOPES

Manufacturer : **XENOSYS CO., LTD.**

2, Venture-ro, 100 beon-gil, Yeonsu-gu, Incheon 406 840, Republic of Korea

Tel. +82 32 875 9811 Fax. +82 32 875 9813

E-mail : sales@xenosys.co.kr / hadlee@hanmir.com URL : http://www.xenosys.co.kr

FEATURES

- Larger field of view & superior depth of focus
- Lighter & Smaller Telescopes
- Clear & Bright vision over the full field
- No Distortion image
- No Fatigue on eye
- IP Frames
 - Super light weight & comfortable aircraft aluminum front
 - Hypoallergenic titanium temple
 - Ion Plating Color

• IPHL Frame



Navy

Brown

Grey(Graphite)

• IPHX Frame



Navy

Brown

Grey(Graphite)

SPECIFICATIONS

Model	Mag.	F.O.V. @40cm	Weight*
Looks2500	2.5x	160 mm @ 40 cm	33.0g
Looks2500L	2.5x	150 mm @ 40 cm	30.0g
Looks2800	2.8x	130 mm @ 40 cm	31.8g
Looks3200	3.2x	86 mm @ 40 cm	49.1g
Looks3500	3.5x	98 mm @ 40 cm	69.1g
Looks5000	5.0x	62 mm @ 40 cm	77.4g
Looks7500	7.5x	41 mm @ 40 cm	83.6g
Flip-up	Magnification : 2.5x Working Distance choice among: 30 cm, 36 cm, 42 cm, 48 cm		

* Weight based on TTL loupe with plano lens with IPAL frame

Flip-up 2.5x



FDA CE CFDA ISO13485:2016 ISO14001:2015

Looks 2500



Looks 2800



Looks 3200



Looks 3500



Looks 5000



NEW Looks 7500



Authorised Distributors : JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072, e-mail : jollytechnocrats@gmail.com, jollytechnocratsservice@gmail.com

website : www.jollytechnocrats.in, Mobile : +91-9810025581, +91-9650595509

- ✓ Maximum Brightness 110,000 lx. (LLMX566D model, Battery powered @40cm)
- ✓ Single LED X566 have 61,000 lx. (Distance @ 40 cm away)
- ✓ Magnetic extension cable: Easy to attach & detach, comfortable to store
- ✓ Comes with light-weight handy carrying case • Battery running time maximum 24 hours (at the lowest level)

New Launch

Additional battery (INCLUDED)



Remote Controller



Carrying Case



L2SX5-S (Single LED)				
Step	Brightness (lx) @WD400mm Illumination (±10%)			Operating Time
	LLMX582	LLMX574	LLMX566	
1 st	12,000	14,000	18,000	24H
2 nd	22,000	25,000	33,000	13H
3 rd	30,000	40,000	44,000	10H
4 th	41,000	47,000	61,000	6h30m
L2SX5 -D(Double LED)				
Step	Brightness (lx) @WD400mm Illumination (±10%)		Operating Time	
	LLMX566D			
1 st	36,000		12h50m	
2 nd	64,000		7h30m	
3 rd	87,000		5h00m	
4 th	110,000		3h10m	

FEATURES

- Super bright Dual LED headlight even its battery powered with ONE control unit
- Compact & light weight control unit, LED module
- Improved & compact size battery pack & Secure assembling system
- LED heat conduction effect: Multi link & attachment is made with aluminum
- LED module weight 12.3 g (LLMX582)
- More soft & durable LED wires

 www.xenosys.co.kr	Manufacturer: XENOSYS CO., LTD.
	2, Venture-ro, 100 beon-gil, Yeosu-gu, Incheon 406-840, Republic of Korea
	Tel: +82 32 875 9811 Fax: +82 32 875 9813
	E-mail: sales@xenosys.co.kr / hadlee@hanmir.com URL: http://www.xenosys.co.kr

Authorised Distributors :

JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072, e-mail : jollytechnocrats@gmail.com, info@jollytechnocrats.in

website : www.jollytechnocrats.in, Mobile : +91-9810025581, +91-9650595509

