

E200ix Trauma Strecher



Specifications

E200ix Trauma Stretcher

This stretcher with entire X-Ray function designed for whole body. It also added 90 degrees sitting posture and another function is that installed the pneumatic cylinder to adjust crook for knees from 0 degree to 40 degree.

- + Integrated storage tray with oxygen bottle trough
- + Trend / Reverse Trend
- + Double side mounted foot control hydraulic pump
- + X Ray trans lucnt
- + Brake pedal control four wheels / steering pedal in two end of bed
- + IV pole inserted hole x 5
- + Protective corner bumpers x4
- + Collapsible siderails
- + Pneumatic assist backrest bends the angle $0^{\circ} \sim 90^{\circ}$
- + Retractable fifth-wheel steering
- + 4" (10cm) mattress or more
- + 8" width wheel
- + Wheel cover
- + Delivery push bar
- + Two stages movable IV rod x1
- + Two stages permanent folding IV Rod x1
- + Foot auxiliary push bar
- + No-crank pneumatic knee gatch $0^{\circ} \sim 40^{\circ}$
- + Upright oxygen holder aperture

