

p27

C500 Operating Chair



Standard Features

Adjustable contoured backrest

- + Dual-density foam seat cushions
- + Foot-controlled, hydraulic, three-spindle base
- + Front-end locking system
- + ABS Thermoformed plastic hood
- + Weight Capacity 300 lbs (136kg)



Seat Cushion Style

	Water drop Height Range
High	70 cm±2
Low	53 cm±2

	Wedge Height Range
High	72 cm±2
Low	55 cm±2

	Saddle Height Range
High	70 cm±2
Low	53 cm±2

PSI Maximus
Multi purpose
Advanced Electro Hydraulic
Surgical table

Variants
PSI Maximus TS
With Top Slide
PSI Maximus OBZ
Bariatry

PSI Legend
Multi purpose
Electro Mechanical
Time tested Surgical Table

Variants
PSI Legend Ortho
Dedicated Ortho Table
PSI Legend TS
With Top Slide
PSI Legend Uro
Dedicated Urology Table
PSI Legend Pedo
Dedicated Pediatric Table
PSI Legend OBG
Labor cum Obstetric
surgical table.
PSI Legend Litho
Lithotripsy table
PSI Legend Eye
Dedicated Eye Surgery table

PSI Virtue
Multi Purpose
Mechanical Surgical Table

Variants
PSI Virtue Ortho
Dedicated Ortho Table

PSI Levitate
Foot operated
Motorized & Hydraulic
Surgeons Chairs

PSI Durus
Heavy Duty
Patient Recovery Trolley
With stretcher

PSI Atlas
Motorized & Static
Operating Room Pendants

PSI Sterilis
Pressure Steam sterilizers

PALAKKAD SURGICAL INDUSTRIES PVT. LTD.

Petronet Road, Wise Park, Kanjikode East
Palakkad, Kerala, India 678 621
Tel: +91 491 2569000, 2566443, 2566551
Fax: +91 491 2569911
Email: palakkadsurgical@gmail.com
www.palakkadsurgical.com

PSI-06/2012

PSI Levitate

Foot operated **Surgeons Chairs**

PSI
PALAKKAD SURGICAL INDUSTRIES PVT. LTD.

A synonym, for quality **surgical equipments** over the decades.

We are overwhelmed by the persistence of our customers for Excellence.

We have always felt that we are in an environment where excellence is not only expected but demanded.

Commitment to excellence

For years we had been relentlessly pursuing excellence and it has now become a habit, rather an attitude for us.

We are prepared, competent, comfortable to meet your demands for excellence.

We Offer You
PSI Levitate

PSI Levitate HD

Hydraulic surgeon's chair

Surgeon while sitting can raise/lower the height Hydraulically by pedaling of foot pump.

Minimum height 23" and maximum 32"

Comfortable for use of microscope while sitting

Height adjustable / detachable arm rests.

Improvised versions of arm rests can be made to order

Adjustable back rest has vertical and horizontal movements.

Full swiveling castors for easy mobility. Castors with brake optional.

Complete stainless steel parts with upholstered seat



PSI Levitate MT

Motorized surgeon's chair

Surgeon while sitting can raise/lower the height electrically by foot switch

Minimum height 20" and maximum 28". Comfortable for use of microscope while sitting

Height adjustable / detachable arm rests.

Improvised versions of arm rests can be made to order

Adjustable back rest has vertical and horizontal movements.

Full swiveling castors for easy mobility. Castors with brake optional.

Complete stainless steel parts with upholstered seat

