



Phillips MX30 Prescription Radiation Glasses - Part # RG-MX30

SKU: RG-MX30-RX

The MX30 frame is the **first protective sports eyewear to receive the AOA seal of acceptance**. The MX30 **Meets ANSI Z87.1+ standards** along with ASTM F803 **impact resistance** standards for several different applications.

It has a **unisex shape** and **temple bars and nose bridges that have soft inside rubberized padding** for extra comfort.

It comes with semi-ridged case and adjustable nylon strap.

Featuring our standard **0.75mm lead equivalency Schott SF6 radiation safety glass lenses**.

Available colors: Gloss Black, Navy Blue, Red and Silver

Models offered:

Model
RGMX30-PLANO/NON PRESCRIPTION
RG MX30-SV PRESCRIPTION
RGMX30-FE BIFOCAL PRESCRIPTION
RG-PROGRESSIVE PRESCRIPTION



123 Lincoln Blvd. • Middlesex, NJ • 08846

Phone: 732-356-1493 • Fax: 732-356-7127

Email: bbrown@phillips-safety.com