



FITOVER RADIATION GLASSES

MODEL RG-33



RG-33-BK (Black)



RG-33-W (White)



RG-33-T (Tortoise)

DESCRIPTION

Phillips Model RG-33 fitover radiation protective glasses offer a comfortable, lightweight, and durable design that fits over prescription eyewear. The frame is constructed of our durable TR-90 nylon and offers .35mm Pb lateral protection. This product offers our standard high quality, distortion-free SF-6 Schott glass radiation-reducing lenses with .75mm lead equivalency. The Model RG-33 glasses will fit over almost any style of glasses and head sizes with comfort.

FEATURES

- *0.75mm Lead Equivalency Schott SF6 Radiation Protective Lenses*
- *Ultimate Fitover Design for Weight & Comfort*
- *0.35mm Pb Lateral Protection*

SPECS

Frame Size: 62-15-145
Weight: 90g
CE 0120 Certified
Rx Frame

