

## Blood Donor Couch DC-100X & DC-60X (Blood Donor Couch (DC-100X), Portable Blood Donor Lounge)

**Description :** DC-100X Blood Donor Couch is a completely automatic enveloping variable tilt chair. DC-100X was specifically designed to make blood withdrawals easier, safe and functional and it them turned to be perfect also for other diagnostic and therapeutic areas. The possibility to taking a semi orthopneic position, which is undoubtedly preferred by patients of blood donors compared with the dorsal position lying on a bed makes long blood transfusion chemotherapeutical of hemodialysis sessions more bearable. Immediate rotation. In case of lipothymia a lightpressure on the control button makes the whole chair immediately tilt along its transversal axis, thus lowering the head and raising the legs, in a moderate Head low & foot high position. As a consequence, it favors a rapid blood flow to the brain districts and an immediate recovery from the lipothymial state. Slowly and gradually the position is then brought back to the initial semi-seated and then upright position, in order to avoid orthostatic hypotension events. Equipped with a practical I.V. rod, with a height-adjustable hook, that can be easily installed on the hard ware. These characteristics make the DC-100X couch the ideal solution for Blood Transfusion Centers, Diagnostic Analysis Labs and Day Hospitals.



### Other Features:

- Provides a comfortable position for the donor
- Variable positioning for either arm with a comfortably wide arm-rest
- Arm rest have swinging out as well as up and down moving facility
- Reclining and upright body positions with a smooth shifting to any position
- Both sides have supporting brackets
- Comfortable working level for the operator
- Lifting capacity – Approx 200 kg.
- Drawers provided for the upkeep of equipment.

### Features:

- Seat, head rest & arm rest can be adjusted to the phlebotomy postures.
- Elector-mechanical control of movements enable head high/low position.
- Suitable height for convenience of phlebotomist.
- Remote Switch for easy repositioning.
- Lockable castors for easy mobility.
- Swing in/out arm rests.
- 100% St. Steel framing & Hardware.
- Washable Upholstery & Punishment proof laminates.

### Specification:

External Dimension	220 cms / Full Length, Seat Length – 60 cms, Back Rest – 92 cms
Loading weight	200 Kg
Voltage Requirements	220-240 V AC
Average Power Input	60 Watts
Arm Rest	Left/right side adjustable with UP/Down Movement 18 cms wide
Scope of Delivery	AC connection cables, Arm Rests ** Musical Device
Movement Control	Push Up / Push Down
Weight	90 kgs. Approx
Packing	One unit (Transport worthy pack)
Time for Position Change	45 Seconds