

Full-Featured Systems for PAD Diagnosis, Satisfies Reimbursement Requirements

Vista AVS™

Faster, Easier Segmental Exams – A Full-Featured ABI/Vascular System

This full-featured vascular assessment system is ideally suited to perform segmental exams. The Vista AVS offers 3 modalities for obtaining systolic pressures and waveforms required for reimbursement and is designed to accommodate both standard and custom protocols. The system includes:

- 3 modalities: PVR, bi-directional Doppler probe and PPG probe
- Automatic ABI, TBI and segmental calculation
- Performs seated ABI on mobility impaired patients
- Software for full-page reporting, exam storage and interface with EMR
- 9 blood pressure cuffs with automatic inflation/deflation
- Waveform display provides easy-to-follow instructions

OPTIONAL FEATURES FOR THE VISTA AVS

5MHz Bi-Directional Probe

For deeper vessel penetration on larger patients, this probe can be selected at time of purchase in addition to the 8MHz bi-directional. The SD5V probe car also be ordered for use with the Vista ABI

Case for Vista AVS/ABI, Vantage ABI

Sold separately, this carrying case conveniently holds the main unit and its accessories for easy transport.

Cuff Selector

This segmental cuff selector is specially designed to make segmental pressure measurements easy to perform. Up to 10 vascular cuffs can be connected to this unit at the same time, for upper and lower extremity or digit pressure measurements.



ITEM NUMBER	DESCRIPTION
L500VA	A Full-Featured ABI/Vascular System ^Δ
SD5V	5 MHz Bi-Directional Probe
K280	Case for Vista AVS/ABI, Vantage ABI
K290	Cuff Selector

The Vista ABI

This Vista ABI system was designed for practices that want to begin performing ABI exams now and also have the option of later adding accessories like the PPG probe and software. This system includes:

- 2 modalities: PVR and bi-directional Doppler probe
- Automated cuff inflation/deflation
- Built-in printer
- Graphic waveform display with step-by-step ABI instructions
- 4 cuffs
- Optional PPG package and software package

The Vista ABI does not accommodate segmental exams.

OPTIONAL FEATURES FOR THE VISTA ABI

PPG Package

A sensitive arterial PPG probe for obtaining both toe and limb pressures. Comes with convenient Velcro to secure placement on the digit. This package includes a PPG probe and a 1.9cm cuff for performing TBI exams on patients with non-compressible arteries.

PC Download Software

An optional software package for the Vista ABI with full-page reporting and exam storage includes download capability for integration with EMR systems (.pdf and .csv file formats – USB Interface to PC).



ITEM NUMBER	DESCRIPTION
L450VA	Vista ABI System ^Δ
L450PG	PPG Package
L450SF	PC Download Software