

# Vascular Assessment Dopplers and Probes That Easily Detect Vessels and Blood Flow

## LifeDop™ 150 Series

These handheld, non-display Dopplers for vascular assessment offer superior sound quality and function with the full line of non-directional LifeDop probes. Rechargeable batteries and audio recording options are available on certain models.

ITEM NUMBER	DESCRIPTION
L150-SD4	LifeDop 150 with 4 MHz Probe <sup>Δ†</sup>
L150-SD5	LifeDop 150 with 5 MHz Probe <sup>Δ†</sup>
L150-SD8	LifeDop 150 with 8 MHz Probe <sup>Δ†</sup>
L150-STN	LifeDop 150 with 8 MHz Sterilizable Probe <sup>Δ†</sup>
L150A-SD4	LifeDop 150 with Audio Recorder & 4 MHz Probe <sup>Δ†</sup>
L150A-SD5	LifeDop 150 with Audio Recorder & 5 MHz Probe <sup>Δ†</sup>
L150A-SD8	LifeDop 150 with Audio Recorder & 8 MHz Probe <sup>Δ†</sup>
L150R-SD4	LifeDop 150 with 4 MHz Probe, Rechargeable <sup>Δ†</sup>
L150R-SD5	LifeDop 150 with 5 MHz Probe, Rechargeable <sup>Δ†</sup>
L150R-SD8	LifeDop 150 with 8 MHz Probe, Rechargeable <sup>Δ†</sup>
L150R-STN	LifeDop 150 with 8 MHz Sterilizable Probe, Rechargeable <sup>Δ†</sup>
K200	Roll Stand with Storage Basket
SINPK1	L150 Doppler with Sterilizable Probe and Stand

## LifeDop™ 250 Series

Offering rich features, these handheld Dopplers have large displays for heart rate and battery status. Rechargeable batteries and audio recording options are available on specific models. They accept all LifeDop non-directional probes.

ITEM NUMBER	DESCRIPTION
L250-SD4	LifeDop 250 Disp. with 4 MHz Probe <sup>Δ†</sup>
L250-SD5	LifeDop 250 Disp. with 5 MHz Probe <sup>Δ†</sup>
L250-SD8	LifeDop 250 Disp. with 8 MHz Probe <sup>Δ†</sup>
L250-STN	LifeDop 250 Disp. with 8 MHz Sterilizable Probe <sup>Δ†</sup>
L250AR-SD4	LifeDop 250 Disp. with Audio Rec., Recharge & 4 MHz Probe
L250AR-SD5	LifeDop 250 Disp. with Audio Rec., Recharge & 5 MHz Probe
L250AR-SD8	LifeDop 250 Disp. with Audio Rec., Recharge & 8 MHz Probe
L250AR-STN	LifeDop 250 Disp. with Audio Rec., Recharge & 8 MHz Sterilizable Probe <sup>Δ†</sup>
L250R-SD4	LifeDop 250 Disp. with Recharge & 4 MHz Probe <sup>Δ†</sup>
L250R-SD5	LifeDop 250 Disp. with Recharge & 5 MHz Probe <sup>Δ†</sup>
L250R-SD8	LifeDop 250 Disp. with Recharge & 8 MHz Probe <sup>Δ†</sup>
L250R-STN	LifeDop 250 Disp. with Recharge & 8 MHz Sterilizable Probe <sup>Δ</sup>
K200	Roll Stand with Storage Basket

## LifeDop™ 350 Series

The LifeDop 350 is a convenient tabletop Doppler with enhanced features for more efficient detection and blood flow monitoring. It has a large heart rate display, rechargeable batteries, and accepts all LifeDop non-directional probes. Enhanced features include:

- Excellent sound/high fidelity speaker
- Improved ergonomics and overall appearance
- Two probe holders for multiple applications
- Louder audio output
- User accessible batteries
- Improved stand and wall mounting

ITEM NUMBER	DESCRIPTION
L350R-SD4	LifeDop Tabletop Doppler with 4 MHz Probe <sup>Δ†</sup>
L350R-SD5	LifeDop Tabletop Doppler with 5 MHz Probe <sup>Δ†</sup>
L350R-SD8	LifeDop Tabletop Doppler with 8 MHz Probe <sup>Δ†</sup>
K310	Wall Mount for Added Convenience and Security
K220	Convenient Roll Stand with Storage Basket
STNPK2	L350R Doppler with Sterilizable Probe and Stand

## LifeDop™ Vascular Assessment Probes

All Summit Doppler LifeDop vascular probes are specifically designed to provide excellent sensitivity for locating vessels and blood flow. Extremely durable, LifeDop probes are interchangeable for use on the L150, L250, and L350 series of Dopplers, with the exception of the sterilizable probes, which have some restrictions. Most models have user replaceable cables.

### 4 MHz Vascular Probe

This probe is optimized with a broad, flat beam for easy location of blood flow.

### 5 MHz Vascular Probe

The medium width beam aids in detecting deep vessels.

### 8 MHz Vascular Probe

The best general-purpose vascular probe. The sensitive, narrow beam allows for easy detection of surface vessel location.

### 8 MHz Sterilizable Probe

Clinically tested for use in sterile fields for the detection of blood flow in both veins and arteries. STERRAD® compatible for use with selected LifeDop Dopplers.

ITEM NUMBER	DESCRIPTION
SD4	4 MHz Vascular Probe <sup>†</sup>
SD5	5 MHz Vascular Probe <sup>†</sup>
SD8	8 MHz Vascular Probe <sup>†</sup>
STN	8 MHz Sterilizable Vascular Probe <sup>†</sup>
C250CC	Obstetrical Vascular Combination – L250 Doppler, 3 MHz and 8 MHz Probe
C250CV	Vascular Combination – L250 Doppler, 5 MHz and 8 MHz Probe <sup>Δ</sup>



SummitDoppler™  
A WALLACH DEVICE



SummitDoppler™  
A WALLACH DEVICE