

DW-BL5 Aluminum Alloy Foldaway Stretcher 折叠担架(铝合金)



四折担架


Technical features:

- The stretcher is made of high-strength aluminum alloy and oxford leather materials. It is characterised by its being light-weight, small-sized, easy-carried, use-safely and easy for sterilization. It's mainly used for hospitals, sports, ambulance and battle fields carrying patients and wounded person.

技术特点:

- 本担架采用高强度铝合金材料和牛津革制成，具有重量轻，体积小，携带方便、使用安全、易消毒，清洗等优点，主要适用医院、公共场所、救护车及部队战地运送伤员。

DW-BL6 Aluminum Alloy Foldaway Stretcher 折叠担架(铝合金)



四折担架

Technical features:

- The stretcher is equipped of outrigger, portable backpack, made of high-strength aluminum alloy and Oxford materials.
- It characterised by its being light-weighted, small-sized and easy-carried, Use-safety and easy for sterilization It's mainly used for field rescue.

技术特点:

- 本担架带支撑腿，带手提背包。采用高强度铝合金焊接而成，台面为牛津革和特斯林。
- 它具有重量轻，体积小，携带方便、使用安全、易消毒清洗等优点，主要适用于野外救护。

Model型号	Product Size 展开尺寸 (L x W x H)	Foldaway 折叠尺寸 (L x W x H)	Net weight 自重	Load Limit 承重	Packing 纸箱包装		
					Quantity 数量	Size 尺寸(L x W x H)	G.W.毛重
DW-BL5	214 x 52 x 12cm	58 x 21 x 16cm	5kg	160kg	1	62 x 28 x 17cm	5.5kg
DW-BL6	202 x 54 x 16cm	54 x 22 x 20cm	8kg	160kg	1	56 x 24 x 22cm	8.5kg

DW-BS1 Stainless Steel Foldaway Stretcher 折叠担架(不锈钢)



Technical features:

- This stretcher is intensified stretcher, it can be folded with flexible handle, outrigger.
- It is made of high-strength stainless steel materials. It characterised by its being small-sized use safely and easy for sterilization and clean.
- It's mainly used for battle field and outdoor carrying patients and wounded person.

技术特点:

- 本担架为加强型担架，可折叠，配有手柄，带支腿。
- 采用高强度不锈钢材料制成，它具有体积小、使用安全、易消毒清洗等优点，主要适用于部队野外救护。

DW-BS2A Stainless Steel Basket Stretcher 篮式担架(不锈钢)

DW-BS2B Steel Basket Stretcher 篮式担架(碳钢)



Technical features:

- The stretcher is made with an eye on a very wide range of special emergency situations, such as in the mountains, air and at sea.
- It is sturdy and flexible to use, its quick reliable fittings enable first-aid person to operate it rapidly and safely. With its special sling equipment, the stretcher is ideal for lifted and transported by helicopter.
- It has an adjustable feet-securing mechanism, two safety belts and a mattress.
- The basket stretcher is made of stainless steel tube and net, it will not rust and potentially cause additional harm to the patient. It is stronger than common basket stretcher.
- It can be used together with the Spine Boards, into two parts, easy to carrying and storage.

技术特点:

- 此担架的构造着眼于急救的特殊性，如崎岖的山区、空中或海上救援。
- 框架坚固耐用，简便可靠的装置让操作人员能够安全快捷的采取急救措施。
- 悬吊能与飞机上挂构连接，实现野外救援。担架配有安全带和垫子。
- 由不锈钢管(碳钢)和网制成，不会对病人造成另外的伤害，比普通担架更牢固，可与板式担架一起使用；可分成两部分，易运输和储存。

Model型号	Product Size 展开尺寸 (L x W x H)	Foldaway 折叠尺寸 (L x W x H)	Net weight 自重	Load Limit 承重	Packing 纸箱包装		
					Quantity 数量	Size 尺寸(L x W x H)	G.W.毛重
DW-BS1	220 x 52 x 12cm	110 x 52 x 15cm	6kg	160kg	1	120 x 54 x 18cm	7kg
DW-BS2A	219 x 69 x 19cm	120 x 69 x 20cm	17kg	380kg	1	125 x 71 x 22cm	25kg
DW-BS2B	219 x 69 x 19cm	120 x 69 x 20cm	17kg	380kg	1	125 x 71 x 22cm	25kg