



DL2 Hybrid M



DL3N



DL4

	DL2 Hybrid M	DL3N	DL4
Polarization	Cross & Non	Cross & Non	Cross & Non
Infection Control (IceCap™)	-	-	Yes
Spectrum Control	-	PigmentBoost	PigmentBoost Plus
Optics	ø25 mm 10x	ø25 mm 10x	ø30 mm 10x
Battery	Rechargeable	Rechargeable	Rechargeable
Camera-compatible	Yes	Yes	Yes
Phone & iPad-compatible	Yes	Yes	Yes
10 mm markings	Yes	Yes	Yes
Included accessories	Silicone sleeve	Desktop charger Leather pouch Eyepiece	USB cable Belt Holster Silicone Sleeve Eyepiece IceCap sample supply

LT-2015B | Specifications subject to change without notice.



dermlite DL4

3Gen®

©2001-2015 by 3Gen, Inc. All rights reserved. Printed in the United States | For information about 3Gen patents & trademarks, please visit www.dermLite.com/patents
 3Gen, Inc. | 31521 Rancho Viejo Road Suite 104 | San Juan Capistrano, CA 92675 | U.S.A. | Tel. +1-949-481-6384 | Fax +1-949-240-7492 | info@3Geninc.com | www.dermLite.com
 iPad and iPhone are reg. trademarks of Apple Inc., registered in the US and other countries.



Magnetically quick & simple
DermLite DL4 can be magnetically attached to a DermLite iPhone 6 case without any additional parts.

- + compact all-aluminum design
- + excellent ergonomics
- + instant-on & auto shutoff
- + large 30 mm 10x lens system
- + charge from any USB port
- + polarized & non-polarized
- + PigmentBoost™ Plus
- + contact & non-contact dermoscopy
- + magnetically attach cameras, smartphones & iPad
- + IceCap™ infection control system
- + four charge indicator levels
- + two brightness levels
- + magnetic faceplate
- + sturdy protective belt holster



IceCap™
Disposable snap-on infection control



Belt Holster
Protects your DL4 when on the go.

The most refined vision in dermoscopy.

The 4th-generation DermLite dermatoscope has been designed from the ground up to be more refined in every detail. An amazingly smooth, ergonomically optimized all-aluminum design that comfortably molds to your hand, DermLite DL4 puts phenomenal versatility, the largest optics of any pocket-size DermLite, and instant on/off control right at your fingertips. Built upon the DL3's proven polarized & non-polarized lighting and solid battery life, DL4 now offers two brightness levels, a larger field of view, enhanced imaging of deeper pigmentation with two levels of color spectrum control (PigmentBoost Plus), and is still more than 25% smaller than the previous generation.



DL4 uses magnetic connectivity to achieve groundbreaking convenience and fast workflows. Not only are both faceplate (with 10 mm reticle) and eyepiece magnetically attached, your smart devices, like iPhone, iPad or Galaxy S series phones, also attach magnetically, enabling you to go between taking clinical and dermoscopic photos quicker than ever before. Along with the growing need for dermoscopy, infection control is becoming an increasingly important aspect of providing care. With its Infection Control and Epidemiology system IceCap™, DL4 has you covered. Disposable snap-on IceCaps attach to the spacer and may greatly reduce the risk of cross-contamination.



A newly developed DL4 retractable glass faceplate design allows you to go between non-contact & contact modalities smoother than ever and lets you precisely focus your image while in skin contact. DermLite DL4's high-powered LEDs and rechargeable battery produce light in an extremely energy-efficient way. And, in order to know when you run out of battery, DL4 now has a more granular four-level charge indicator. If you need to recharge, you can do so from any USB port using the standard Micro USB to USB cable included with each DL4. In its elegance and compactness, DermLite DL4 is simply the most capable dermatoscope available.