


## 5. Models



---

### CryoIQ DERM PLUS

Liquid  
Freezing



CIQ-DP-L

Contact  
Freezing



CIQ-DP-C1



CIQ-DP-C3



CIQ-DP-C5

---

### CryoIQ PRO

Liquid  
Freezing



CIQ-P-L

Contact  
Freezing



CIQ-P-C1



CIQ-P-C3



CIQ-P-C5



CIQ-P-C7

## 5.1 Options for Model PRO



Long 5° tip for Liquid "Spray" Freezing

Glass tip, length 130 mm  
Treatment area  $\varnothing$ 7-18 mm

L13/5



Angled tip 45° for Liquid "Spray" Freezing

Glass tip, length 80 mm  
Treatment area  $\varnothing$ 7-18 mm

A45



D1



D2



D3

Tips for Liquid "Spray" Freezing

Glass tip, length 22 mm  
D1= Treatment area  $\varnothing$ 1-6 mm  
D2= Treatment area  $\varnothing$ 7-18 mm  
D3= Treatment area  $\varnothing$ 9-22 mm



C1



C3



C5



C7

Tips for Contact Freezing

Contact surface  $\varnothing$ 1, 3, 5, 7 mm



Cartridge 16 /25 gram with valve

Cartridge designed to work with the  
CryoIQ® | DERM, PRO

Packed in box with 1 cartridge  
Packed in box with 4 cartridges