

DERMATOME ST-5

Surgical Motor Systems

Battery Operated

Power, Precision and Innovation...



■ Guides for Graft Carrier ■ Ratchet Handle for Easy Operation ■ Graft Thickness Adjustment (0-1,5 mm)



CARRIER

Skin Graft Mesher



M11-500G
Skin Graft Mesher
205 x 165 x 120 mm / 4000 g



M11-600 Skin graft carrier > 1,5:1 (Box of 20)
M11-605 Skin graft carrier > 3:1 (Box of 20)

- Single cutter cuts skin graft into geometrical shaped mesh & enlarges skin surface
- Graft expansion predetermined by graft carrying plates, eliminates the need to change the cutter during operation
- Continuous feed design prevents graft carrying plate from being loaded incorrectly and allows multiple plates to be used without interruption
- Graft thickness adjustment feature allows expansion on thin and thick grafts
- Easy availability of replacement parts

