

01 ENT Treatment Unit (IU-3000)



Product Feature

POWER SOURCE	220V AC, 50/60Hz	STANDARD COMPONENT
POWER CONSUMPTION	1.5kW	Spray (3ea)
COMPRESSOR MOTOR	0.9~2.0kg / cm ²	Suction (2ea)
SUCTION MOTOR	Max. Vacuum 680mmHG	Ventilation
	Max. pressure 7bar	Post pole (max 2ea, Straight/ Curved)
ANTI-FOG DEVICE	330W	Medicine bottle (6ea)
DEVICE CONTROL SWITCH	Water resistant	Medicine bottle tray (option, 14ea)
	Sound effect	Internal Instrument tray (Heating)
OPTION	Monitor with arm set	External Instrument tray
	Side can	heater
	Sub table	pen light
	Doctor Chair	Outlet (6ea)
	Patient Chair	100mm Can (2ea)
DIMENSION	W 1096 x D 737 x H 890 (1447)	83mm Can (12ea)
WEIGHT	110kg	53mm Can (2ea)

