

DELIVERY TABLE

DT 100

PHOENIX
Phoenix Medical Systems (P) Ltd.



Comfort for the mother
Convenience for the care-provider

The Phoenix DT 100 is a highly flexible delivery table. It has been designed for use in labour rooms and maternity wards. It is used in hospitals & clinics for obstetric, gynaecological, antenatal and colposcopic procedures.

With the DT 100, the mother can adopt any convenient birthing position.

The squat position



The lithotomy position on calf support



The lithotomy position on foot rest



The DT 100 is available in a choice of different pleasing colours to match the interior of your labour room.



DISCLAIMER: Actual product colour may slightly differ

The doctor, midwife and nurse appreciate the versatility of the DT 100.



The supine position



Motorized adjustable height



STANDARD



U-CUT SEAT



FOOT REST



CALF SUPPORT



CATCH BASIN



FIRM HAND GRIP



SIDE RAILS



FOOT CONTROL



PAPER ROLL HOLDER



ROTATABLE IV POLE



CASTOR WITH BRAKE



ARM REST



SOFT CUSHION BASE

OPTIONAL

- ▶ Hand - operated controls in lieu of foot - operated controls
- ▶ Battery backup supporting upto 50 movements
- ▶ Anaesthesia screen holder
- ▶ Emergency CPR positioning (-10°)
- ▶ Squat bar
- ▶ Reusable bed cover
- ▶ Baby receiving tray



CDSCO Registration No.: 000387 / 2020

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Chennai - 600044, India.
Toll Free : 1 800 425 3503, Tel.: 91-44-2478 1705 / 2227,
✉ India: sales@pmsind.com, International: exports@pmsind.com
🌐 : phoenixmedicals.com