

# Bella LED

A powerful lamp on a very robust flexible arm appreciated by specialists who wish to perform small interventions in their practice, and who need skin and tissue tone rendered accurately. The ergonomic handle enables precise and intuitive positioning.

The optional intensity switch (50/100 %) allows the practitioner to view the skin under different light levels.

Ideal for dermatology, and minor or cosmetic surgery

Medical device class 1

- Light intensity **1680 Lumens\***
- Colour of temperature **4500°K\* (neutral white)**
- Colour rendering index **CRI 95/100**
- Electric consumption **17 Watts**

Lux*	Beam 10°	Distance	Illuminated diameter
210000		30 cm	10 cm
120000		50 cm	13 cm
33000		100 cm	20 cm

\* Values accurate within +/- 10 %



Bella  
contactless switch

- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2:41
- Power cable length: 3 meters
- Suggested colors: white
- Head dimensions: 16 x 13 cm

- Fittings: 1, 2, 3, 4, 6, 7, 8, 9, 11, 12, 13

CODES		
	On/off Switch	Contactless switch
Flexible 70 cm	LED17650	LED17650-HSC
Flexible 70 cm Double intensity option	LED17650DI	LED17650DI-HSC



# Hepta LED

Professional examination lamp with a timeless design that combines power, solidity, and versatility for heavy use.

With a tough flexible neck and an ergonomic handle, it can be positioned as required and fulfils all the needs of medical lighting.

The optional dimmer switch allows the brightness to be adjusted for user comfort and optimal use.

Ideal for all types of use in general medicine, dermatology, gynaecology, urology, podiatry, etc.

Medical device class 1

- Light intensity **980 Lumens\***
- Colour of temperature **4500°K\* (neutral white)**
- Colour rendering index **CRI 95/100**
- Electric consumption **7 Watts**

Lux*	Beam 10°	Distance	Illuminated diameter
130000		30 cm	10 cm
60000		50 cm	12 cm
18000		100 cm	20 cm

\* Values accurate within +/- 10 %



Hepta 65  
Dimmer switch

- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2:41
- Power cable length: 3 meters
- Suggested colors: white and black
- Head dimensions: 11 x 10 cm

- Fittings: 1, 2, 3, 4, 6, 7, 8, 9, 10, 11, 12, 13

CODES			
	On/off Switch	Dimmer switch	Contactless switch
Flexible 65 cm	LED07651	LED07651V	LED07651-HSC
Flexible 114 cm	LED071001	LED071001V	LED071001-HSC

