

JMM

03

**JMM 03****NITROCARE ÜÇ MOTORLU JİNEKOLOJİK MUAYENE MASASI**  
NITROCARE 3 MOTORS GYNAECOLOGICAL EXAMINATION COUCH**JMM 003**

Jinekolojik ve jübolojik muayene ve ufak cerrahi işlemler için uygundur. Üç motorlu sırt, yükseklik ve translemburg hareketlerini el kumandası aracılığıyla gerçekleştirmektedir.

Muayene masası baş/sart – oturum – ayak olmak üzere üç ayrı bölüme ayrılmaktadır.

İntegral malzemenin imalı, ayarlanabilir, isenildiği zaman tabele çıkarılabilir, ayılbabilir destek mekanizması mevcuttur.

Yükseklik ayarlaması, ağırlıkla destekli, ayarlanabilir, kolay temizlenebilir, pelillerin sünger üzeri aynı dertli saphmalı örtürme – yama yüzeyi mevcuttur.

Masa low-bed özelliğine sahiptir ve yıkama/ısıtma hareketi esnasında pozisyonunda herhangi bir değişime uğramaz.

**Sırt Kısmı Dik Konumdayken:**

- Maks: 1800 mm
- Min: 1290 mm

**Sırt Kısmı Yatık Konumdayken:**

- Maks: 1010 mm
- Min: 500 mm

Patlamaz metal malzemenin imal 9'lü çeyre çeyre tasarılan, sahip olduğu teleskopik ray yardımıyla ileri geri kesonelli hareket ettirilebilen ağırlık saphmalı hareket mekanizması.

**JMM 003**

JMM 003 examination table is designed for gynaecological and urological examinations and minor operations.

It has 3 motors enabling back, height and translemburg adjustments via remote control. It consists of 3 sections: head/back, seating and feet.

Additionally, it comes with adjustable and removable foot/calf support which is made of integral material.

Lying surface is made of antibacterial, water resistance and easy to clean polyurethane foam covered by faux leather.

Examination table has low-bed feature and keeps its position during height adjustments

**Back at vertical position**

- Max: 1800 mm
- Min: 1290 mm

**Back at horizontal position**

- Max: 1010 mm
- Min: 500 mm

It includes a 9 litre capacity stainless steel waste disposal bucket which sits on a rail enabling backward and forward movement.

DMO KODU 355.685.086

LOW BED



SIRT HAREKETİ



YÜKSEKLİK



TRENDELEMBURG

Ağırlık Toplama Kuvveti  
Stainless Steel, Removable Bstein