



- 1-) 360° Dönebilen Tekerler
- 2-) Ayarlanabilir Baldır Destek Aparatı
- 3-) Mekanik CPR
- 4-) Açı Ölçer
- 5-) Korkuluk Açma Kapama
- 6-) Çarpma Koruyucu Tamponlar
- 7-) El Kumandası
- 8-) Serum Askılığı

- 1-) 360° Movable Castors
- 2-) Adjustable Leg Stirrups
- 3-) Mechanical CPR
- 4-) Protractor
- 5-) Rails Control (Open-Close)
- 6-) Plastic Protection Bumpers
- 7-) Hand Panel
- 8-) Height Adjustable IV Pole

Merkezi Fren Sistemi  
Central Brake SystemPaslanmaz Ahtık Toplama Küveti  
Stainless Steel, Removable BasinÇarpma Koruyucu Tamponlar  
Plastic Protection BumpersAçı Ölçer  
Protractor

## TEKNİK ÖZELLİKLER

JMM02 Jinekoloji yatağı muayene üroloji operasyon ve hastanın operasyon sonrası dinlenebileceği konforda üretilmiştir.

Sırt ve yükseklik hareketlerini motor yardımı ile yapabilmektedir.

Oturma ve yatma yüzeyi Compact malzemeden imal.

Manuel kontrollü CPR

Manuel Trendelenburg

Uluslararası standartlara sahip yan korkuluk koruması.

360° dönebilen merkezi kilitli tekerlekler (125mm)

Yüksekliği ve pozisyonu ayarlanabilen ayak destekleri

Paslanmaz çelik leğen

Poliüretan malzemeden imal kilitlenebilen yan korkuluklar ve çıkarılabilen karyola başlıkları.

Plastik koruyucu çarpma tamponları

İdrar torbası askılığı

Kemer Kancası

Aksesuar Yuvaları

Yüksekliği ayarlanabilir serum askısı

Elektrostatik toz boya

32 dansilite şilte

## TECHNICAL SPECIFICATIONS

Excellent Use For Examinations, Urological And Simple Operations And Post Operations  
Backrest And Height Adjustments By Motors

HPL Compact Lying Platform

Mechanical CPR

Mechanical Trendelenburg

Side Rails Protection For Patients As Per International Standarts

125 mm Castors With Central Brake System

Height And Position Adjustable Leg Stirrups

Stainless Steel, Removable Basin

Lockable Side Rails And Detachable Heading Produced From PU ABS

Plastic Protection Bumpers

Urine Bag Hanger

Belt Hooks

Accessorie Slots

Height Adjustable IV Stand

Electrostatic Powder Coated Metal Frame And Components

32 Density Mattress

