

LEEP PRECISION™ Integrated System

Loop Electrosurgical Excision Procedures



Fully integrated and mobile, the LEEP Precision System combines three components – Cart, Generator and Smoke Evacuator.



LEEP PRECISION Integrated System

NUMBER	DESCRIPTION	
LP-10-120	LEEP PRECISION Integrated System, 120V	CA*
LP-10-220	LEEP PRECISION Integrated System, 220V International	CE

LEEP PRECISION Generator

NUMBER	DESCRIPTION	
LP-20-120	LEEP PRECISION Generator, 120V	CA*
LP-20-220	LEEP PRECISION Generator, 220V International	CE

LEEP PRECISION Smoke Evacuator

NUMBER	DESCRIPTION	
LP-30-120	LEEP PRECISION Smoke Evacuator, 120V, Stand-Alone	CA*
LP-30-220	LEEP PRECISION Smoke Evacuator, 220V, Stand-Alone International	CE

LEEP PRECISION Cart

NUMBER	DESCRIPTION	
LP-40-001	LEEP PRECISION Cart	CA* CE

LEEP PRECISION Generator Disposables

NUMBER	DESCRIPTION	
LP-50-101	Hand Switch-Operated Push-Button Cautery Pencils (10/Box)	CA*
LP-50-102	Foot Pedal-Operated Cautery Pencils (10/Box)	CA*
LP-50-103	Foot Pedal Adapter (1/Bag)	CA*
LP-50-201	Dual Dispersive Electrodes with Connector (10/Box)	CA*
LP-50-202	Dual Dispersive Electrodes without Connector (Connect to Reusable Cable) (10/Box)	CA*
LP-50-203	Reusable Cable (1/Bag)	CA*



LEEP PRECISION Smoke Evacuator Disposables

NUMBER	DESCRIPTION	
52560	Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag)	CA*
6082	Smoke Evacuator ULPA Filter (1/Box)	CA*
6083	Smoke Evacuator Reducer, Single-Patient Use (10/Bag)	CA*
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)	CA*

LEEP Electrodes

(Shown Actual Size)

Our broad selection of 21 different offerings includes loops, squares, needle and ball tips specifically designed for loop electrosurgical excision procedures.



Large-Radius Loop Electrodes ≥ 2 cm Wide CA* CE

NUMBER	DESCRIPTION
R2510	2.5 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2020	2.0 cm W x 2.0 cm D • 12 cm shaft (5/Box)
R2015	2.0 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R2010	2.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2008	2.0 cm W x 0.8 cm D • 12 cm shaft (5/Box)

Medium-Radius Loop Electrodes <2 cm Wide CA* CE

NUMBER	DESCRIPTION
R1515	1.5 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R1507	1.5 cm W x 0.7 cm D • 12 cm shaft (5/Box)
R1505	1.5 cm W x 0.5 cm D • 12 cm shaft (5/Box)
R1010	1.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R1007	1.0 cm W x 0.7 cm D • 12 cm shaft (5/Box)

Square Electrodes CA* CE

NUMBER	DESCRIPTION
S1010	1.0 cm W x 1.0 cm D 12 cm shaft (5/Box)
S1007	1.0 cm W x 0.7 cm D 5.5 cm shaft 15° bend (5/Box)
S1004	1.0 cm W x 0.4 cm D 5.5 cm shaft 15° bend (5/Box)
V1004	1.0 cm W x 0.4 cm D 8 cm shaft (5/Box)

Needle Electrodes CA* CE

NUMBER	DESCRIPTION
N1200	0.8 mm D x 16 mm 12 cm shaft (5/Box)
N0550	0.8 mm D x 16 mm 5.5 cm shaft (5/Box)
N2512	0.25 mm D x 8 mm 12 cm shaft (5/Box)
N2555	0.25 mm D x 5 mm 5.5 cm shaft (5/Box)

Ball Electrodes CA* CE

NUMBER	DESCRIPTION
B0512	5 mm D 12 cm shaft (5/Box)
B0555	5 mm D 5.5 cm shaft (5/Box)
B0312	3 mm D 12 cm shaft (5/Box)


Fischer Cone Biopsy Excisor®

For Cone Biopsy Procedures

The Fischer Cone Biopsy Excisor (FCBE) comes in a variety of sizes and shapes.

Fischer Cone Biopsy Excisor CA* CE

NUMBER	DESCRIPTION
900-150	Fischer Cone Biopsy Excisor, Small (5/Box)
900-151	Fischer Cone Biopsy Excisor, Medium (5/Box)
900-152	Fischer Cone Biopsy Excisor, Large (5/Box)
900-154	Fischer Cone Biopsy Excisor, Medium Extended (5/Box)
900-155	Fischer Cone Biopsy Excisor, Large Shallow (5/Box)
900-156	Fischer Cone Biopsy Excisor, Mixed Set (S, M, ME, L & LS) 1 each (5/Box)
900-157	Fischer Cone Biopsy Excisor, Small Wide-Angle (5/Box)
900-158	Fischer Cone Biopsy Excisor, Small Wide-Angle Extended (5/Box)



LEEP Procedure Essentials



LEEP RediKit®

For sale in the USA only.

Each single-use LEEP RediKit contains the following:

- 1 Vial Xylocaine (Lidocaine HCl 1%) 10 ml with Epinephrine*
- 1 Potocky Needle® 27 G (for administration of cervical block)*
- 1 Lugol's (Strong Iodine Solution, USP) 8 ml, Single Dose
- 1 AstrinGyn® (Ferric Subsulfate), 8 grams, Single Dose
- 6 Applicator Swabs, 8 inch
- 20 G Needle
- 12 cc Luer Lock Syringe*
- 1 Pantyliner
- 1 Alcohol Wipe*
- 1 Underpad

*Sterile

NUMBER	DESCRIPTION
6061	LEEP RediKit, Single-use, disposable. Contact Customer Service for ordering information.



Needle Extension

NUMBER	DESCRIPTION
433-132	Needle Extension • 10 cm (1/Box)
433-133	Needle Extension • 13 cm (1/Box)
433-134	Needle Extension • 15 cm (1/Box)



Lugol's

Lugol's iodine solution. For sale in the USA only.

NUMBER	DESCRIPTION
6064	Lugol's. Contact Customer Service for ordering information. (12/Box)
6064-4	Lugol's. Contact Customer Service for ordering information. (48/Case)



AstrinGyn®

Monsel's solution. For sale in the USA only.

NUMBER	DESCRIPTION
6065	AstrinGyn (12/Box)
6065-4	AstrinGyn – 4 Boxes of 12 (48/Case)



Potocky Needle®

27-gauge Potocky Needle. The Potocky Needle attaches to a "dental style" cartridge syringe, while the Potocky Needle with Luer Lock is used with a standard medical syringe.

NUMBER	DESCRIPTION
6066	Potocky Needles (6/Box)
6066-4	Potocky Needles – 4 Boxes of 6 (24/Case)
6166	Potocky Needles with Luer Lock (6/Box)
6166-4	Potocky Needles with Luer Lock – 4 Boxes of 6 (24/Case)

LEEP System Accessories



Foot Pedal, LEEP System 1000

Featuring pneumatic, air-operated control, the foot pedal does not use electric current.

NUMBER	DESCRIPTION
6032	Foot Pedal for LEEP System 1000




Return Pad, Single-Patient Use CA*

Non-sterile patient return pads.

NUMBER	DESCRIPTION
6050P1	LEEP Return Pad, Single-Patient Use (10/Box)

LEEP Accessories



Compatible Disposable Accessories

See your local sales representative for more information.

NUMBER	DESCRIPTION
395-210	Electrode Tip Cleanser, Disposable (24/Box)
395-302	Conductive Electrode Pads, Disposable for ESU (25/Box)
395-305	Conductive Electrode Pads, Dual Dispersive, Disposable (25/Box) Plus II Gen.
395-310	Hand-Switch Pencil and Cord, Disposable (25/Box)
72308	Dual Section Grounding Pads (10/Box)
72102	Filter/Hose Set for Smoke Evacuator (Pack of 3)
72311	Patient Return Cable
22192	LEEP 6000 Return Pad
22191	Foot Pedal



Tubing Set, Single-Patient Use CA*

Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuator) tube. This set completes the connection between the reducer and the specula.


NUMBER	DESCRIPTION
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)



Handpiece, Single-Patient Use CA*

The sterile, single-patient use handpiece works with any of our disposable electrodes. The wand is operated with a footswitch.

NUMBER	DESCRIPTION
6040	Handpiece, Single-Patient Use (10/Box)



Disposable Smoke Evacuator (DSE) Tubing

Non-sterile single-use, DSE tubing is designed to work specifically with the LEEP rose-coated specula. The 9-inch DSE tubing is flexible and noncombustible and contains no harmful PVC.

Also available as a set with 6 feet of tubing. See Tubing Set, Single-Patient Use, product item number 6084.

NUMBER	DESCRIPTION
6150	Disposable Smoke Evacuator Tubing – 9" straw (50/Bag)
34371	LEEP System 1000 Manual (No charge)

Patient Return Pad Adapter LEEP System 1000



This reusable replacement adapter connects the single-use Patient Return Pad to the LEEP System 1000.

NUMBER	DESCRIPTION
23635	Patient Return Pad Adapter for LEEP System 1000 (for use with 6050P1 pads)



Reducer CA*

Disposable, single-patient use reducer connects the Pre-Filter (52560) with the Patient Tubing Set (6084) on the Smoke Evacuator.

NUMBER	DESCRIPTION
6083	Smoke Evacuator Reducer (10/Box)



Pre-Filter CA*

Replacement filters for the CooperSurgical Smoke Evacuator. Used to eliminate smoke odor.

NUMBER	DESCRIPTION
52560	Smoke Evacuator Pre-Filter (10/Box)



ULPA Filter CA*

Used with the Smoke Evacuator.

NUMBER	DESCRIPTION
6082	Smoke Evacuator ULPA Filter (1/Box)