


GYNAEC & GENERAL APPLICATIONS

MANUAL HAND VACCUM PUMP

For minor work with required safety negative pressure achieved. Suction pressure can be controlled or seen visually and the user can stop or proceed with the function accordingly since it is controlled by ourselves (Hand pump). Very low cost. No electricity required and hence can be moved to any place or for camp work. Silicon Vaccum Cup size available small / medium / large.

**Any make silicon cup can be attached.*



This **SAFETY FILTER** should be used before procedure to avoid liquid flow inside the pump valve.