



### Consistent Pressure

VacuSense™ Technology only in the Ameda Platinum™ pump automatically adjusts to single and double pumping and to atmospheric pressure at different locations for a consistent pumping experience.



### Fully Adjustable

CustomControl™ allows complete control of suction and speed to select best levels for milk flow and comfort.



### Sizes To Fit

CustomFit Flange System™ provides a range of breast flange sizes to help make pumping more comfortable and efficient.

## Ameda Breast Pumps



### Una™ Single Electric Breast Pump

Perfect for mothers who prefer single pumping.

- Independent customizable speed and suction controls.
- Proven Airlock Protection reduces risk of contamination for safe milk.



The simple and effective choice for mothers returning to work or who want a quick, easy and effective breast pumping experience.

### Lactaline™ Double Electric Breast Pump

- Separate suction and speed dials allow for CustomControl™ Dual adjustability lets the mother select the best levels for her milk flow and comfort, creating a multi-phase pumping experience
- Airlock Protection™ prevents airflow between expressed milk and pump tubing while pumping
- Convenient power source options: AC adapter, AA batteries or optional car adapter



#### Includes:

- Lactaline Electric Breast Pump
- Double Hygiene Pumpset
- AC adapter



### Manual Breast Pumps

Compact, ergonomic one-hand manual pumps earned the High Design-Innovation Quality Award from Design Zentrum Nordrhein-Westfalen Germany.

Easy-to-use, allowing you to achieve multi-phase pumping by varying handle squeezes for more milk flow.



## DISCOVER SAFE MILK & MORE MILK

### WITH THE AMEDA DIFFERENCE

Ameda's clinically proven breast pumps and breastfeeding accessories have provided nursing mothers with an effective, comfortable breastfeeding experience since 1942.



#### Safe Milk

Proven Airlock Protection™ with an FDA-cleared barrier that helps protect milk against bacteria and viruses by preventing milk getting into the tubing.



#### More Milk

ComfortFlow™ Technology creates a smooth, more consistent sensation so pumping feels comfortable, which means more milk. Shown effective for establishing and maintaining adequate milk production in hospitals.

#### Authorised Distributors:

### JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of: OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 80, Nehru Place, New Delhi-110019, (India)

Phone: +91-11-46066071, Telefax: +91-11-46066072, e-mail: jollytechnocrats@gmail.com, jollytechnocratservice@gmail.com

website: www.jollytechnocrats.in, Mobile: +91-9810520055



## Ameda Milk Collection Systems

Ameda has a clinically proven protective barrier that helps keep milk and moisture out of the pump tubing and motor while pumping.

You can enjoy the benefits of using the same kit with all Ameda Electric Breast Pumps.

**The Ameda Difference**



**Ameda Barrier**



All Ameda pumps are designed as a closed system with an FDA-cleared barrier to help keep milk safe.

**Non-Ameda**



Mold can be a dangerous consequence of not having an effective and proven barrier.

**MOLD**

<sup>1</sup> Jones, E. et al. A randomized controlled trial to compare methods of milk expression after preterm delivery. Arch Dis Child Fetal Neonatal Ed 2007;90: F95.



### Ameda Cool 'N' Carry Milk Storage Bag

Ameda Cool 'N' Carry™ Milk Storage bag with 3 freezer elements For discreet cooling and transport of breast milk, use one Cool 'N' Carry Tote for daycare and another at work, or have two to alternate when transporting breast milk for a hospitalized baby. Insulated zipper nylon zipper case can be easily cleaned by wiping with a damp cloth. Integrated handle for easy carrying. Discreet appearance.

Holds up to six 4 oz./120 mL bottles and three cooling elements. Includes:

- (1) Insulated black nylon carry bag
  - (3) Cooling elements
  - (1) Milk storage guidelines card
  - (1) Instructions for use
- Bottles are BPA free and DEHP free

## CustomFit™ Flanges

Every mother needs a good breast flange fit for greater comfort and better milk flow.

Research has found that nearly one third of mothers who pump report greater comfort and more milk with a larger-than-standard breast flange size.<sup>1</sup>

Ameda CustomFit Flanges are available in seven sizes and have been designed to work with all Ameda Milk Collection Systems.



NO BPA

## Purely Yours Ultra



- Proven Airlock Protection™ barrier helps protect milk while keeping pump tubing dry and clean. No need to clean tubing.
- CustomControl™ empowers moms to more closely mimic their baby's nursing rhythm. Moms have complete control of their cycle, both suction and speed, to select the best levels for their milk flow and comfort.

#### Purely yours Ultra

Purely Yours Ultra (is the same as PY Lactaline), Dottie Tote(bag), European power adapter, Dual Hygienikit® Milk collection system, 2 valves, custom fit breast flanges (2) 30.5mm & 28.5mm reducing inserts, cool & carry Tote (1) insulated bag, 6(4oz) bottles with lock-tight lids, (3) cooling elements.



**Distribution in other countries:** For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit [www.ameda.com](http://www.ameda.com).