

UV COMB

Simple and effective unit with narrow band UVB Lamp in the treatment of psoriasis, vitiligo, specially on the scalp and localized parts of the body. The lightweight unit ergonomically developed, easy to handle. The unit can also be used without the comb for partial irradiation of small areas. Short exposure enables time saving application. other special features. The extensive device contains all the parameters to successfully treat, Psoriasis, Vitiligo, Atopic Dermatitis, Circumscribed Scleroderma and other skin disorders. It also offers highest level of safety.



Specialty of UV comb:

Use of NBUVB lamp.

Parabolic aluminium reflector e n s u r e s maximum output.

On / Off switch.

Separate comb attachment for s c a l p treatment.

Technical data:

Type: UV Comb (UV-B-311 nm)
Power supply: 230 V. AC / 50 Hz
Lamp: PL-S 9W/01 (Philips)
Dimensions (L x W x H): 270 x 50 x 50 mm
Radiation outlet: 160 x 25 mm
Operating temperature: 5 - 45 °C

