




LEEP PRECISION™ Integrated System

Loop Electrosurgical Excision Procedures



Fully integrated and mobile, the LEEP Precision System combines three components – Cart, Generator and Smoke Evacuator.

LEEP PRECISION Integrated System

NUMBER	DESCRIPTION	CA#
LP-10-120	LEEP PRECISION Integrated System, 120V	CA#
LP-10-220	LEEP PRECISION Integrated System, 220V International	CE

LEEP PRECISION Generator

NUMBER	DESCRIPTION	CA#
LP-20-120	LEEP PRECISION Generator, 120V	CA#
LP-20-220	LEEP PRECISION Generator, 220V International	CE

LEEP PRECISION Smoke Evacuator

NUMBER	DESCRIPTION	CA#
LP-30-120	LEEP PRECISION Smoke Evacuator, 120V, Stand-Alone	CA#
LP-30-220	LEEP PRECISION Smoke Evacuator, 220V, Stand-Alone International	CE

LEEP PRECISION Cart

NUMBER	DESCRIPTION	CA#
LP-40-001	LEEP PRECISION Cart	CA# CE


LEEP PRECISION Generator Disposables

NUMBER	DESCRIPTION	CA#
LP-50-101	Hand Switch-Operated Push-Button Cautery Pencils (10/Box)	CA#
LP-50-102	Foot Pedal-Operated Cautery Pencils (10/Box)	CA#
LP-50-103	Foot Pedal Adapter (1/Bag)	CA#
LP-50-201	Dual Dispersive Electrodes with Connector (10/Box)	CA#
LP-50-202	Dual Dispersive Electrodes without Connector (Connect to Reusable Cable) (10/Box)	CA#
LP-50-203	Reusable Cable (1/Bag)	CA#



LEEP PRECISION Smoke Evacuator Disposables

NUMBER	DESCRIPTION	CA#
52560	Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag)	CA#
6082	Smoke Evacuator ULPA Filter (1/Box)	CA#
6083	Smoke Evacuator Reducer, Single-Patient Use (10/Bag)	CA#
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)	CA#




LEEP

LEEP Electrodes

(Shown Actual Size)

Our broad selection of 21 different offerings includes loops, squares, needles and ball tips specifically designed for loop electrocautery excision procedures.



Large-Radius Loop Electrodes ≥ 2 cm Wide CA* €€

NUMBER	DESCRIPTION
R2510	2.5 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2020	2.0 cm W x 2.0 cm D • 12 cm shaft (5/Box)
R2015	2.0 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R2010	2.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2008	2.0 cm W x 0.8 cm D • 12 cm shaft (5/Box)

Medium-Radius Loop Electrodes <2 cm Wide CA* €€

NUMBER	DESCRIPTION
R1515	1.5 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R1507	1.5 cm W x 0.7 cm D • 12 cm shaft (5/Box)
R1505	1.5 cm W x 0.5 cm D • 12 cm shaft (5/Box)
R1010	1.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R1007	1.0 cm W x 0.7 cm D • 12 cm shaft (5/Box)

Square Electrodes CA* €€

NUMBER	DESCRIPTION
S1010	1.0 cm W x 1.0 cm D 12 cm shaft (5/Box)
S1007	1.0 cm W x 0.7 cm D 5.5 cm shaft (5/Box)
S1004	1.0 cm W x 0.4 cm D 5.5 cm shaft 15° bend (5/Box)
V1004	1.0 cm W x 0.4 cm D 8 cm shaft (5/Box)

Needle Electrodes CA* €€

NUMBER	DESCRIPTION
N1200	0.8 mm D x 1.6 mm 12 cm shaft (5/Box)
N0550	0.8 mm D x 1.6 mm 5.5 cm shaft (5/Box)
N2512	0.25 mm D x 8 mm 12 cm shaft (5/Box)
N2555	0.25 mm D x 5 mm 5.5 cm shaft (5/Box)

Ball Electrodes CA* €€

NUMBER	DESCRIPTION
B0512	5 mm D 12 cm shaft (5/Box)
B0555	5 mm D 5.5 cm shaft (5/Box)
B0312	3 mm D 12 cm shaft (5/Box)


Fischer Cone Biopsy Excisor®

For Cone Biopsy Procedures

The Fischer Cone Biopsy Excisor (FCBE) comes in a variety of sizes and shapes.

Fischer Cone Biopsy Excisor **CA** **CE**

NUMBER	DESCRIPTION
900-150	Fischer Cone Biopsy Excisor, Small (5/Box)
900-151	Fischer Cone Biopsy Excisor, Medium (5/Box)
900-152	Fischer Cone Biopsy Excisor, Large (5/Box)
900-154	Fischer Cone Biopsy Excisor, Medium Extended (5/Box)
900-155	Fischer Cone Biopsy Excisor, Large Shallow (5/Box)
900-156	Fischer Cone Biopsy Excisor, Mixed Set (S, M, ME, L & LS) 1 each (5/Box)
900-157	Fischer Cone Biopsy Excisor, Small Wide-Angle (5/Box)
900-158	Fischer Cone Biopsy Excisor, Small Wide-Angle Extended (5/Box)



LEEP Procedure Essentials



LEEP RediKit®

For sale in the USA only.

Each single-use LEEP RediKit contains the following:

- 1 Vial Xylocaine (Lidocaine HCl 1%) 10 ml with Epinephrine*
- 1 Potocky Needle* 27 G (for administration of cervical block)*
- 1 Lugol's (Strong Iodine Solution, USP) 8 ml, Single Dose
- 1 AstrinGyn® (Ferric Sub sulfate), 8 grams, Single Dose
- 6 Applicator Swabs, 8 inch
- 20 G Needle
- 12 cc Luer Lock Syringe*
- 1 Partyliner
- 1 Alcohol Wipe*
- 1 Underpad

*Sterile

NUMBER	DESCRIPTION
6061	LEEP RediKit, Single-use, disposable. Contact Customer Service for ordering information.



Needle Extension

NUMBER	DESCRIPTION
433-132	Needle Extension • 10 cm (1/Box)
433-133	Needle Extension • 13 cm (1/Box)
433-134	Needle Extension • 15 cm (1/Box)



Lugol's

Lugol's iodine solution. For sale in the USA only.

NUMBER	DESCRIPTION
6064	Lugol's. Contact Customer Service for ordering information. (12/Box)
6064-4	Lugol's. Contact Customer Service for ordering information. (48/Case)



AstrinGyn®

Monsef's solution. For sale in the USA only.

NUMBER	DESCRIPTION
6065	AstrinGyn (12/Box)
6065-4	AstrinGyn – 4 Boxes of 12 (48/Case)



Potocky Needle®

27-gauge Potocky Needle. The Potocky Needle attaches to a "dental style" cartridge syringe, while the Potocky Needle with Luer Lock is used with a standard medical syringe.

NUMBER	DESCRIPTION
6066	Potocky Needles (6/Box)
6066-4	Potocky Needles – 4 Boxes of 6 (24/Case)
6166	Potocky Needles with Luer Lock (6/Box)
6166-4	Potocky Needles with Luer Lock – 4 Boxes of 6 (24/Case)

LEEP System Accessories



Foot Pedal, LEEP System 1000

Featuring pneumatic, air-operated control, the foot pedal does not use electric current.

NUMBER	DESCRIPTION
6032	Foot Pedal for LEEP System 1000



Handpiece, Single-Patient Use

The sterile, single-patient use handpiece works with any of our disposable electrodes. The wand is operated with a footswitch.

NUMBER	DESCRIPTION
6040	Handpiece, Single-Patient Use (10/Box)



Patient Return Pad Adapter LEEP System 1000

This reusable replacement adapter connects the single-use Patient Return Pad to the LEEP System 1000.

NUMBER	DESCRIPTION
23635	Patient Return Pad Adapter for LEEP System 1000 (for use with 6050P1 pads)



Return Pad, Single-Patient Use

Non-sterile patient return pads.

NUMBER	DESCRIPTION
6050P1	LEEP Return Pad, Single-Patient Use (10/Box)



Disposable Smoke Evacuator (DSE) Tubing

Non-sterile single-use, DSE tubing is designed to work specifically with the LEEP rose-coated specula. The 9-inch DSE tubing is flexible and noncombustible and contains no harmful PVC.

Also available as a set with 6 feet of tubing. See Tubing Set, Single-Patient Use, product item number 6084.

NUMBER	DESCRIPTION
6150	Disposable Smoke Evacuator Tubing – 9' straw (50/Bag)
34371	LEEP System 1000 Manual (No charge)

LEEP Accessories



Compatible Disposable Accessories

See your local sales representative for more information.

NUMBER	DESCRIPTION
395-210	Electrode Tip Cleanser, Disposable (24/Box)
395-302	Conductive Electrode Pads, Disposable for ESU (25/Box)
395-305	Conductive Electrode Pads, Dual Dispersive, Disposable (25/Box) Plus II Gen.
395-310	Hand-Switch Pencil and Cord, Disposable (25/Box)
72308	Dual Section Grounding Pads (10/Box)
72102	Filter/Hose Set for Smoke Evacuator (Pack of 3)
72311	Patient Return Cable
22192	LEEP 6000 Return Pad
22191	Foot Pedal

Tubing Set, Single-Patient Use

Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuator) tube. This set completes the connection between the reducer and the specula.

NUMBER	DESCRIPTION
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)



Reducer

Disposable, single-patient use reducer connects the Pre-Filter (52560) with the Patient Tubing Set (6084) on the Smoke Evacuator.

NUMBER	DESCRIPTION
6083	Smoke Evacuator Reducer (10/Box)



Pre-Filter

Replacement filters for the CooperSurgical Smoke Evacuator. Used to eliminate smoke odor.

NUMBER	DESCRIPTION
52560	Smoke Evacuator Pre-Filter (10/Box)



ULPA Filter

Used with the Smoke Evacuator.

NUMBER	DESCRIPTION
6082	Smoke Evacuator ULPA Filter (1/Box)



Euro-Med™ LEEP Instruments

These nonconductive, fully autoclavable LEEP instruments are offered in a variety of styles to facilitate cervical tissue excision and removal.




Pederson LEEP Specula

NUMBER	DESCRIPTION	SIZE MM
F223	Small	W13 L81
F221	Medium	W22 L114
F220	Large	W25 L120



Vu-More™ LEEP Specula

NUMBER	DESCRIPTION	SIZE MM
F217	Medium	W33 L113
F215	Large	W35 L127



Vu-Max™ LEEP Specula

NUMBER	DESCRIPTION	SIZE MM
F222	Pederson Medium	W22 L114
F224	Pederson Large	W25 L120
F218	Graves Medium	W33 L113
F219	Graves Large	W35 L127



Mini Vu-More™ LEEP Specula

NUMBER	SIZE MM
F250	W34 L93

(W=Distal End) (L=Penetration Length)



Open-Sided LEEP Specula

NUMBER	SIZE MM
F210	W35 L107

Prima Series® Vaginal LEEP Specula **CA***

NUMBER	DESCRIPTION	SIZE MM
PSV-1L	Medium	W33 L113
PSV-2L	Large	W35 L125



Graves LEEP Specula

NUMBER	DESCRIPTION	SIZE MM
F205	Medium	W33 L113
F200	Large	W35 L127



Snowman™ LEEP Specula

NUMBER	DESCRIPTION	SIZE MM
F107	Graves Blade	W35 L152
F117	Pederson Blade	W25 L152
F127	One each/Graves and Pederson	




Pederson LEEP Specula (with Integrated Smoke Tube)

NUMBER	DESCRIPTION	SIZE MM
61521	Small	W16 L84
61518	Medium	W22 L114
61519	Large	W26 L121



Graves LEEP Specula (with Integrated Smoke Tube)

NUMBER	DESCRIPTION	SIZE MM
61513	Small	W22 L87
61514	Medium	W33 L113
61515	Large	W35 L125





Euro-Med® LEEP Instruments

Tru-View™ LEEP Lateral Wall Retractor

NUMBER	DESCRIPTION	SIZE MM
F317	Lateral Wall Retractor	W32 L190




Jackson LEEP Retractor

NUMBER	DESCRIPTION	SIZE MM
F235	Medium	W40 L89
F234	Large	W40 L108



Cer-View™ LEEP Lateral Wall Retractor

NUMBER	DESCRIPTION	SIZE MM
F410	Medium	W22 L142
F420	Large	W22 L155





LEEP Tenacula, Hooks and Rulers



LEEP Ring Forceps

NUMBER	INSTRUMENT LENGTH
F920	220 mm



LEEP Tissue Forceps

NUMBER	INSTRUMENT LENGTH
F700	210 mm



LEEP Champion Forceps

NUMBER	INSTRUMENT LENGTH
F925	220 mm



LEEP Three-Prong Hook

NUMBER	INSTRUMENT LENGTH
F803	260 mm



LEEP Iris Hook

NUMBER	INSTRUMENT LENGTH
F800	263 mm



LEEP Hook

NUMBER	INSTRUMENT LENGTH
F810	Straight 250 mm
F812	Left-Angled 250 mm
F813	Right-Angled 250 mm



LEEP Depth Ruler

NUMBER	INSTRUMENT LENGTH
F600	240 mm



LEEP Width Ruler

NUMBER	INSTRUMENT LENGTH
F610	210 mm



LEEP Emmett Tenaculum

NUMBER	INSTRUMENT LENGTH
F805	250 mm



LEEP Two-Prong Hook

NUMBER	INSTRUMENT LENGTH
F802	260 mm



LEEP Single-Tooth Tenaculum

NUMBER	INSTRUMENT LENGTH
F900	240 mm

LEEP Endocervical Specula



Large

NUMBER	TIP SIZE	LENGTH
F500	5 mm	273 mm



Medium

NUMBER	TIP SIZE	LENGTH
F510	4 mm	273 mm





LEEP

Electrosurgical System Accessories for Wallach's Quantum 2000*



Hand Switch—Operated Pencils
Sterile, single-use only. For use with Quantum 2000* and other compatible electrosurgical generators. Convenient CUT and COAG push buttons.

NUMBER	DESCRIPTION
909089	Hand Switch—Operated Pencils, Sterile (25 per Box)



Foot Pedal—Operated Pencils

Sterile, single-use only. For use with Quantum 2000* and other comparable electrosurgical generators. Reusable adapter required. Use with any standard disposable electrodes.


NUMBER	DESCRIPTION
909028	Foot Pedal—Operated Pencils, Sterile (25 per Box)



Disposable Split Grounding Pads (with Connector)

For use with Quantum 2000*.

NUMBER	DESCRIPTION
909078	Disposable Split Grounding Pads with Connector (10 per Box)



EndoBlock Needle (for Intra-Cervical Block)

Sterile, single-use only, reinforced with an extended shaft. Threaded hub needle for dental style cartridge syringe. Luer lock for use with medical syringe.

NUMBER	DESCRIPTION
920017	EndoBlock Needle Luer Lock (6 per Box)
920017-1	EndoBlock Needle Luer Lock (60 per Box)
920021	EndoBlock Needle Threaded Hub (6 per Box)



Endocervical Block Kit (for Intra-Cervical Block)

Single-use only. Reinforced long 18 gauge EndoBlock needle for drawing anesthetic. Disposable syringe with finger grip.

NUMBER	DESCRIPTION
920009	Endocervical Block Kit (10 Kits per Box)



Disposable Smoke Evacuator Speculum Tubing (Non-Sterile)

NUMBER	DESCRIPTION
920001	Disposable Speculum Tubing (50 per Box)



Disposable Smoke Evacuator Tubing with Attached Speculum Tubing

Sterile, single-use six-foot-long tubing connects to the smoke evacuator filter assembly.

NUMBER	DESCRIPTION
920002	Disposable Smoke Evacuator Tubing with Speculum Tubing (10 per Box)

