

Login

Username

Customer no.

Password



»Orthopaedics Catalogue »Motor Systems and Hand Drills »Plaster saws

## PLASTER SAWS

Page: 1 2 3

### Eccentric disc

[Art.-No. 40.50.74]



**Description:** Eccentric disc to limit depth of cut for saw blades for plaster saws Oscillo IV, Oscillo V and TIGER


**Characteristics:** for circular saw blades No. 40.50.66 and 40.50.69

**Diameter:** Ø 65 mm

add to the basket

### Circular Saw Blade

[Art.-No. 40.50.68]



**Description:** Circular Saw Blade for plaster saws Oscillo IV, Oscillo V and TIGER

**Characteristics:** for plaster dressing

**Diameter:** Ø 50 mm

add to the basket

### Segment Saw Blade

[Art.-No. 40.50.72]




**Description:** Segment Saw Blade for plaster saws Oscillo IV, Oscillo V and TIGER

**Characteristics:** for plaster dressing

add to the basket

**Depth Saw Blade**

[Art.-No. 40.50.73]



**Description:** Depth Saw Blade for plaster saws Oscillo IV, Oscillo V and TIGER

**Characteristics:** for plaster dressing

add to the basket

**Wrench**

[Art.-No. 40.50.71]



**Description:** Wrench for plaster saws Oscillo IV, Oscillo V and TIGER

**Characteristics:** for exchange of saw blades, pair

add to the basket

**Electric Plaster Saw, complete**

[Art.-No. 40.60.60]



**Author:** OSCILLO Compact

**Description:** Electric Plaster Saw (40.60.61), complete with accessories, 1 saw blade Ø 50 mm (40.62.50), 1 saw blade Ø 65 mm (40.60.65), 1 Depth gauge (40.60.64), in plastic case with power supply

**Characteristics:** **Technical data:** 110 / 240 Volt, 50 / 60 Hz, 100 W

add to the basket

**Electric Plaster Saw, complete**

[Art.-No. 40.60.61]



**Author:** OSCILLO Compact

**Description:** Electric Plaster Saw with power supply

**Characteristics:** **Technical data:** 110 / 240 Volt, 50 / 60 Hz, 100 W

add to the basket

**Circular Saw Blade**

[Art.-No. 40.62.50]



**Description:** Circular Saw Blade for electric oscillating plaster saw OSCILLO Compact

**Characteristics:** for synthetic dressing

**Diameter:** Ø 50 mm

add to the basket

**Circular Saw Blade**

[Art.-No. 40.60.65]



**Description:** Circular Saw Blade for electric oscillating plaster saw OSCILLO Compact

**Characteristics:** for plaster dressing

**Diameter:** Ø 65 mm

add to the basket

**Depth gauge**

[Art.-No. 40.60.64]



**Description:** Depth gauge for electric oscillating plaster saw OSCILLO Compact

add to the basket

Page: 1 2 3