

Haier

Haier Pharmacy Refrigerator Offering Accessibility



HYC-68A

Pharmacy Refrigerator

Reliability

- High efficiency specialty compressor with known field reliability
- Permanently lubricated cooling fan for safety and longevity
- Forced-air cooling with optimized air distribution system design to achieve maximum temperature uniformity and stability
- Optimized refrigeration system design for more effective cooling and speedy recovery

Key Features

- Microprocessor control, LED digital display and temperature adjustment with an increment of 1°C
- Designed for ambient temperature 10-32°C with humidity below 60%

Safety

- Equipped with a complete temperature alarm system featuring audible buzzer and visual flashing light, Power failure
- Capable of alerting failures due to high and low temperature, sensor error, door ajar, power loss
- Safety lock to prevent unauthorized access



HYC-68A



Model	HYC-68A	
Climate Class	N	
Power Supply(V/Hz)	220-240/50	115/60
Power(W)	84	84
Electrical Current(A)	0.65	0.65
Capacity(L/Cu.Ft)	68/2.4	
Net/Gross Weight(approx)	38/40(kg) 83.8/88.2(lbs)	
Interior Dimensions(W*D*H)	415*385*505(mm) 16.3*15.2*19.9(in)	
Exterior Dimensions(W*D*H)	495*580*660(mm) 19.5*22.8*26.0(in)	
Packing Dimensions(W*D*H)	522*620*690(mm) 20.6*24.4*27.2(in)	
Container load (20'/40'/40'H)	108/228/228	
Test Hole	N/A	
Shelves	2	

Haier Medical and Laboratory Products Co., Ltd

Room 403D, Brand Building, Haier industrial park, No. 1 Haier Road, Qingdao, 266101, P.R. China

Tel: +86-0532-88936011 / 7169

Website: www.haiermedical.com

